

ROCKING PLATFORM

Model No : RP - 11

- Capacity : 24 Vaccutainer tube of 4 ml capacity
 or
 2 Tray of Western Blot
- Motion : Tilting - 15 Degree from horizontal
 on both side (Total 30 Degree)
- Main input : 220V / 50Hz Or 110V / 60 Hz
- Speed : 15 rocks / minute
- Platform size: 172 mm X 151 mm
- Weight : 2.2 Kg. Net
- Overall Size : 348 (L) x 198 (W) x 103 (H) mm



Note:- We provide both types of rubber pads with device. Rubber pads are interchangeable.

ROTARY MIXER

Model No : RM - 11

- Mains Input : 230V/50Hz or 110V/60Hz A.C. Supply
- Body : ABS + METAL Body
- Motion : Tilting - fix degree at 38
- Speed : 15 RPM
- Weight : 1.2 Kg.
- Size of Device : 188 (L) x 221 (W) x 260 (H)mm
- Operation : Quiet - No noise
- Standard configuration : 4ml Vaccutainer (EDTA) Tubes - 12 nos
 And ESR TUBE (8mm) - 4 nos.
 (Total-16 nos)
- Special configuration : 16 Tube of either type or mix as required.

